



Brassica transformation




100 % EtOH 2 mins
15 % Sodium hypochlorite 15 mins
3 x SDW




MS, 3% Sucrose, 0.8% phytagar
pH 5.6, plus vitamins

O/N cold
Transfer to C/R for 4 days




Dip petiole
In *Agro* suspension



Co-cultivate

MS, 3% Sucrose,
0.8% phytagar, **2mg/l BAP**
pH 5.6, plus vitamins





72 hrs



x10

Selection

Gamborg's B5
media



glasshouse